

SAMPLE PREPARATION and SOFTWARE ACCESS:

- Obtain Spec. grade CH_2Cl_2 in the hood closest to the instrument. Prepare a concentrated solution of your analyte and apply to a salt plate. Salt plates can be found in a dessicator near the instrument.
- Access the FT-IR software by double clicking on the **Horizon HB** icon on the desktop.
- Enter the following information:

User Name: **FTIR**

Password: **AJM328A.**

- If asked, enter **chemistry** as the department name.

CHECKING SPECTROMETER COMMUNICATION:

- Click on **Acquisition** in the **Instrument** tab of the toolbar.
- Click on the **Connection** button in the Instrument panel.
- Click on **Live** and monitor the instrument signal strength.
- In the **Acquisition Parameter** window, set the **Resolution** according to application needs and adjust the **Detector Gain** value to adjust **Signal Strength** between 25% and 90%.




Figure 1. Live Spectrum Screen.

- Click on **Stop** to exit live mode.

SAVING THE DATA:

- Set the location to save data by clicking on the ... button:




Figure 2. Spectrum Acquisition Screen.

- Enter the **Sample Name** and make sure the **Path** and **Project** folder are defined as follows:
Path: **C:\Users\ABBBomem\Documents\horizon MB Data**
Project: **YOUR FOLDER'S NAME**




Figure 3. Sample Information Editor.

- Check the Data to Project box and select the desired folder.

SPECTRUM ACQUISITION:

- Make sure the sample compartment is empty, and adjust the telescopic purge tubes as close to one another.
- The background will be collected in the single beam mode. The FT-IR will collect the sample spectrum as a ratio, subtracting out the background. In order to collect a background, set the acquisition mode to **Single Beam** and the desired number of **Scans** in the **Acquisition Parameter** window and click on the **Reference** button.
- In order to collect a spectrum, place the IR plate in the sample compartment and set the acquisition mode to **Absorbance** or **Transmittance** click on the **Measure** button.

ANALYSIS:

Baseline Correction:

- Select **Baseline Correction** from the **Mathematics** tab in the toolbar.
- Select the desired baseline correction by adjusting the red square boxes on either end of the correction line. Click on **Calculate** in the mathematics window.
- Close the baseline correction window.




Figure 5. Baseline Correction.

Setting Trace Limits:

- Click on the axis to be modified and apply new limits. These will be reflected in the active window and print outs.

Labeling Peaks:

- Select **Peak Picking** from the **Mathematics** tab in the toolbar. A peak table will appear at underneath the spectrum.

