

INITIALIZING THE INSTRUMENT:

- Turn on the computer. When prompted type the password **cary100**.
- Turn on the instrument power switch located on its left front corner. The lamp should be allowed to warm up for at least 15 minutes before running the instrument.
- Start the software by double-clicking on the **SCAN** icon on the desktop.
- Once the instrument is initialized, the screen should have the **Start** and **Stop** buttons at the top center. If **Connect** appears in place of **Start**, press **Connect** in order bring the instrument online.




Figure 1. Scan Toolbar.

SETTING UP THE EXPERIMENT:

- Press the **Setup** button a dialog box with several tabs will appear.
- On the **Cary** tab, enter in your wavelength units and range (the instrument scans from *high to low* wavelengths), **y range** units and values (the program will autoscale, so these numbers aren't crucial), **average time** (like an integration time, longer values give higher signals), and **data interval** (how closely spaced are your data points, smaller intervals give larger data files). Select **Cycle** if you want to run multiple scans of the same sample.
- On the **Options** tab, enter in the **SBW** (bandwidth, which is usually good to leave at 2 nm unless you're looking for fine structure in your spectrum) and **lamp crossover** wavelength (leave at 350 nm unless expecting important features in the spectrum at that wavelength, and then only move it no more than 10 nm one way or the other).
- On the **Baseline** tab, select **Baseline Correction** if doing a scan with a dilute sample, as the baseline and lamp switching will have a negative effect on the spectrum.

- On the **Reports** tab, enter **Operator Name**. To automatically label peaks on the spectrum, select **Maximum Peak** or **All Peaks** options. The peak threshold limit and labeling options are altered by pressing the **Peak Information** button.
- On the **Auto Storage** tab, select whether to save the scans before or after a run. Selecting after the run (not the default choice) will avoid saving bad data.




Figure 2. Set up Dialog Box.

DATA COLLECTION:

- If running a baseline, insert a blank solvent cuvette in both the reference (rear) and sample (front) cell holders. Press the **Zero** button, then the **Baseline** button to run the baseline. Upon completion of the scan, press the autoscale **Y** button (⏏) to see the baseline. A red baseline should appear in the number box at the upper left corner of the screen, indicating a valid baseline is in memory.
- Place your sample in the sample cell holder and press the **Start** button. You will be asked for a sample name (other than the file name), then the spectrum will start running. Again, **autoscale** the spectrum if it goes offscale. Once the run is completed, a save file box will be presented to allow naming the file. Any peak information will be printed in the report box (the bottom half underneath the plot).
- Text can be added to the plot by pressing the **A** button. When running multiple scans in one session, pressing the traces button (the farthest left button in the toolbar) to select individual data sets can be selected and made visible in the current plot. Please note that all data plots (including any baselines) are still accessible through the 'Graph' menu, and all will be printed out when printing the data. In order to omit any data set, right click on the graph with that data and select **Remove Graph**.

SAMPLE SETTINGS:

