

LAB PROCEDURE:

Green Chemistry in the Undergraduate Organic Laboratory: Microwave-Assisted Synthesis of a Natural Insecticide on Basic Montmorillonite K10 Clay

Matthew R. Dintzner*, Paul R. Wucka and Thomas W. Lyons


Department of Chemistry, DePaul University, Chicago, IL 60614, USA;
mdintzne@depaul.edu

Special Chemical Needs: Montmorillonite K10 clay, saturated aqueous Potassium Carbonate, Sesamol, 3-Methyl-2-butenal, Ethyl acetate, 2N Sodium Hydroxide, saturated aqueous Sodium Chloride, Magnesium Sulfate.

Safety: Safety goggles, gloves, and a laboratory coat or apron should be worn at all times. All experiments must be performed in a fume hood. All chemicals should be considered hazardous and direct physical contact with them must be avoided. If exposed, immediately flush skin with water for at least 10 min, while contaminated clothing is removed. Do not ingest or taste any chemicals under any circumstances.

Student Instructions

Prelaboratory: Search the literature and find the references related to the synthesis of methylenedioxyprococene. Read the references and understand the mechanism pertaining to the synthetic route shown below:



Experimental Procedure:

Synthesis of methylenedioxyprecocene (1)



Using a balance, weigh out 500 mg of Montmorillonite K10 clay. Transfer the clay to a 50 mL Erlenmeyer flask. Using a graduated cylinder measure out 10 mL of saturated aqueous K_2CO_3 and add it to the flask containing the clay. Add a magnetic stir bar to the aqueous clay mixture and stir for approximately 10 min at room temperature. Meanwhile set up a vacuum filtration apparatus (side-arm Erlenmeyer flask, filter adapter, Hirsch funnel with filter paper; connect the side arm of the Erlenmeyer to an aspirator via vacuum tubing). Remove the magnetic stir bar, swirl the aqueous basic clay mixture and collect the clay by vacuum filtration. Allow the clay to dry in the Hirsch funnel for several minutes, then transfer it to a clean watch glass, spread it out with a spatula and allow it to dry in an oven at 110 °C for a period of 1 h. Carefully remove the

hot, basic clay from the oven and allow it to cool to room temperature on the watch glass in a desiccator. While the clay is cooling, weigh out the sesamol (138 mg, 1 mmol) and add it to a 20-mL disposable glass scintillation vial (Wheaton). Using a micropipet draw up the 3-methyl-2-butenal (106 μ L, 1.1 mmol) and add it to the vial with the sesamol. Gently tap the vial to allow the sesamol to dissolve. To the resulting dark solution add all of the cool, dry clay and carefully distribute the mixture evenly using a metal spatula. Place the vial in the center of the carousel platform of a household grade microwave oven, cover the vial with a watch glass or other flat piece of glass, and microwave for a period of 8 minutes (several samples may be microwaved at once). Let the reaction mixture cool to room temperature. Then add 5 mL methylene chloride and gently swirl the vial. Separate the clay from the organic solution by vacuum filtration using a Hirsch funnel. Rinse the clay with two successive 5-mL portions of methylene chloride. Transfer the filtrate to a 250-mL separatory funnel and wash successively with 2 N NaOH (2 x 10 mL), and saturated aqueous NaCl (1 x 10 mL). Dry the organic phase over MgSO_4 . Remove the drying agent by vacuum filtration. Transfer the filtrate to a tarred flask and remove the methylene chloride by gently blowing a stream of air over the crude product until mass is constant. Weigh the flask containing the crude methylenedioxyprococene product (a dark brown, viscous oil) and calculate a percent yield. Analyze the product by TLC, GC-MS, IR, or NMR spectroscopy. Assess the purity of your product and interpret all data.

Exercises

1. List four advantages to using clay as a catalyst for organic reactions.
2. What is meant by the term “green chemistry”?
3. What is the purpose of using a microwave oven for this reaction?
4. Show a mechanism for the base-catalyzed reaction of sesamol with 3-methyl-2-butenal to give MDP.

File : C:\HPCHEM\1\DATA\HCSI177B.D
Operator : MRD
Acquired : 28 Oct 104 4:27 pm using AcqMethod DEFAULT
Instrument : 5971 - In
Sample Name: HCS-I-177-B
Misc Info : AFTER WORKUP W/NaOH
Vial Number: 1



HCS-I-177-B



Current Data Parameters
NAME HCS-I-177-B
EXPNO 1
PROCNO 1


F2 - Acquisition Parameters

Date_ 20041028
Time 18.25
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zg30
TD 65536
SOLVENT CDCl3
NS 16
DS 2
SWH 6172.839 Hz
FIDRES 0.094190 Hz
AQ 5.3084660 sec
RG 406.4
LW 81.000 usec
DE 6.00 usec
TE 0.0 K
D1 1.00000000 sec
MCREST 0.00000000 sec
MCMRK 0.01500000 sec

===== CHANNEL f1 =====
NUC1 1H
P1 12.00 usec
PL1 -1.00 dB
SFO1 300.1318534 MHz

F2 - Processing parameters

SI 32768
SF 300.1300000 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 1.00



HCS-I-177-B

148.22
147.59
141.39
128.18
122.31
114.25
105.63
100.89
99.07

77.43
77.21
77.01
76.59
76.02

27.44



Current Data Parameters
NAME HCS-I-177-B
EXENO 2
PROCNO 1

F2 - Acquisition Parameters
Date_ 20041029

Time_ 14.28
INSTRUM spect
PROBHD 5 mm PABBO BB-
PULPROG zgpg30
TD 65536
SOLVENT CDCl3
NS 1024
DS 4
SWH 17985.611 Hz
FIDRES 0.27439 Hz
AQ 1.8219508 sec
RG 16384
DW 27.800 usec
DE 6.00 usec
TE 0.0 K
D1 2.00000000 sec
d11 0.03000000 sec
MCREST 0.00000000 sec
MCMRK 0.01500000 sec

===== CHANNEL f1 =====
NUC1 13C
P1 8.80 usec
PL1 -2.00 dB
SFO1 75.4752953 MHz

===== CHANNEL f2 =====
CPDPRG2 waltz16
NUC2 1H
PCPD2 90.00 usec
PL2 -1.00 dB
PL12 17.00 dB
SFO2 300.1312005 MHz

F2 - Processing Parameters
SI 32768
SF 75.4677490 MHz
WDW EM
SSB 0
TB 1.00 Hz
GB 0
PC 1.00

